



# THONHAUSER

## CLEANING & VERIFICATION OF DISPENSING SYSTEMS WITH TM DESANA MAX

1. PRE-RINSE CLEANING CONTAINER AND LINE



2. OBSERVE SAFETY ASPECTS!



3. PREPARE SOLUTION AND PULL BLANK



4. FILL THE LINES AND SOAK FOR 3 MIN.



5. CHECK COLOR - REPEAT UNTIL MATCHES BLANK



6. FLUSH OFF WITH TAP WATER, RE-CONNECT KEGS



7. HYGIENE OF COUPLERS AND FAUCETS!




ALWAYS FOLLOW INSTRUCTIONS ON PACKAGE!

24-7 POISONING EMERGENCY HOTLINE: 1-866-706-32-62

[desana-usa.com](http://desana-usa.com)

GRADE YOUR LINES:

Try to get as close to the blank (purple color) as possible.



## How much TM DESANA MAX do I need ?

1 packet TM DESANA MAX = 1.5 gal mixed solution

amount of whole TM DESANA MAX packages  
 volume of mixed TM DESANA MAX solution in gallons, pints or liquid ounces

length of tubing in feet	1 FT	5 FT	10 FT	20 FT	50 FT	100 FT	200 FT	300 FT	400 FT	500 FT	600 FT	700 FT	800 FT	900 FT	1000 FT
3/16"	1/6 oz	3/4 oz	1 3/4 oz	1/4 pt	1/2 pt	1 pt	2 pt	1/2 gal	1/2 gal	3/4 gal	3/4 gal	1 gal	1 gal	1 1/4 gal (1)	1 1/4 gal (1)
1/4"	1/3 oz	1 3/4 oz	1/4 pt	1/2 pt	1 pt	2 pt	1/2 gal	3/4 gal	1 gal	1 1/4 gal (1)	1 1/2 gal (1)	1 3/4 gal (2)	2 gal (2)	2 1/4 gal (2)	2 1/2 gal (2)
5/16"	1/2 oz (1/2)	2 1/2 oz	1/4 pt	3/4 pt	1 1/2 pt	3 1/4 pt	3/4 gal	1 1/4 gal (1)	1 1/2 gal (1)	2 gal (2)	2 1/4 gal (2)	2 3/4 gal (2)	3 1/4 gal (2)	3 1/2 gal (3)	4 gal (3)
3/8"	3/4 oz (1/2)	1/4 pt	1/2 pt	1 pt	2 1/4 pt	1/2 gal	1 1/4 gal (1)	1 3/4 gal (2)	2 1/4 gal (2)	3 gal (2)	3 1/2 gal (3)	4 gal (3)	4 3/4 gal (4)	5 1/4 gal (4)	5 3/4 gal (4)
1/2"	1 1/3 oz	1/2 pt	3/4 pt	1/4 gal	1/2 gal	1 gal	2 gal (2)	3 gal (2)	4 1/4 gal (3)	5 1/4 gal (4)	6 1/4 gal (5)	7 1/4 gal (5)	8 1/4 gal (6)	9 1/4 gal (7)	10 1/2 gal (7)